

# RondelO™ IV-set adapter

|           |                          |    |
|-----------|--------------------------|----|
| <b>EN</b> | Instructions for Use     | 2  |
| <b>DE</b> | Gebrauchsanweisung       | 7  |
| <b>ES</b> | Instrucciones de uso     | 12 |
| <b>FR</b> | Mode d'emploi            | 17 |
| <b>IT</b> | Istruzioni per l'uso     | 22 |
| <b>NL</b> | Gebruiksaanwijzing       | 27 |
| <b>PT</b> | Instruções de utilização | 32 |
| <b>FI</b> | Käyttöohjeet             | 37 |
| <b>SK</b> | Návod na použitie        | 42 |
| <b>IL</b> | הוראות לשימוש            | 47 |
| <b>EL</b> | Οδηγίες χρήσης           | 52 |
| <b>HU</b> | Használati útmutató      | 57 |

## 1. General Information

**Intended purpose:** The RondelO IV-set adapter is intended for the infusion/administration of various commonly used medical fluids.

**Indications:** Patient medical condition in need of infusion.

**Contraindications:** Once an intravenous medication is delivered, it cannot be retrieved. When giving IV medications, there is very little opportunity to stop an injection if an adverse reaction or error occurs. IV medications, if given too quickly or incorrectly, can cause significant harm.




**Patient target groups:** General patient population

**Intended users and environment:** Various hospital departments such as intensive care, emergency rooms and diagnostic units. The use of the device should be limited to the specialist trained to perform the procedures for which this device is intended.

**Medical personnel responsibility for patient safety:** In order to competently manage a patient, the medical personnel are required to:

- respond appropriately (adjust treatments as necessary)
- observe for complications and troubleshoot as needed

**Symbols:** The following symbols are used for important information in the Instructions for Use, on the packaging and on the product.

|   |  |   |   |
|---|--|---|---|
|    | Caution                                |    | Consult Instructions for Use  |
|    | Medical device                         |    | DEHP free   |
|    | Manufacturer                           |    | Does not contain, and there is no presence of, natural rubber latex |
|    | Date and country of manufacture        |    | Do not re-use   |
|    | CE mark of compliance                  |    | Single sterile barrier system. Sterilized using ethylene oxide      |
|    | Use-by date                            |    | Do not re-sterilize   |
|   | Batch code                             |   | Do not use if package is damaged                                    |
|  | Catalogue number                       |  | Keep away from sunlight   |
|  | Unique device identifier               |  | Keep dry  |
|  | Non-pyrogenic                          |  | Handling mark on package – Fragile                                  |
|  | Drops per millilitre                   |  | Handling mark on package - Recyclable material                      |
|  | Liquid filter with 15-micron pore size |  | Handling mark on package – This side up                             |
|  | One-way valve                          |   |   |

## Abbreviations

DC - Drip Chamber

FLL - Female Luer Lock

IV-set - Intravenous infusion set

NFV - Needle Free Valve

MLL - Male Luer Lock

## 2. Warnings and Precautions



- The device is designed and intended for single use only. Do not autoclave, reshape, re-sterilize or reuse. Exposure to high heat may cause changes to the functionality of the device.
- Do not re-sterilize. Improper reuse of the medical device can lead to the potential spreading of life-threatening infections.
- Sets are to be changed every 72 hours.
- Make sure the primary infusion is suitable as a flushing fluid for all drugs that will be used.
- Verify the device material is suitable for all drugs that will be used. This device contains PVC.
- Use aseptic technique and follow hospital standard operating procedures when handling the product.
- Respond appropriately and adjust treatments, as necessary.
- Observe for complications and troubleshoot as needed.
- Visually inspect the disinfected surface for possible cracks due to exposure to alcohol.
- Serious adverse events should be reported to Sky Medical a.s. and the applicable Competent Authority.
- This device is not intended for blood or blood products or enteral solutions.

### RondelO positions

Turn the handle to select the intended position. Each position is verified by a haptic click.

- **Drug positions** – Selectable Drug positions correspond to Luer lock connectors of the RondelO device. When one Drug position is selected, the RondelO is open at the corresponding Luer lock connector, and is closed at all other Luer lock connectors.
- **Flushing positions** – Selectable Flushing positions are located in between the Drug positions to allow primary infusion to flush solely and rinse the RondelO internally.

Note: Primary infusion is available in all RondelO positions. Concurrent flushing from Primary infusion is available by default in all Drug positions. Availability of Primary infusion is controlled by the Clamp above RondelO or with head height differential.


 It is highly recommended to consult the instructions prior to use.

Fig 1. Priming of the drip chamber

## 3. RondelO IV-set adapter

**DEVICE DESCRIPTION:** The RondelO IV-set adapter may consist of various components: (1) Vented Spike with protective covering, (2) Pinch clamp, (3) One-way valve, (4) RondelO with 2, 4, or 6, Drug positions at a corresponding Luer lock connector with either a (5) Protective cap, or (6) Needle free valve (NFV), (7) RondelO Handle, (8) Spike port, Tubing, Transparent or UV-protected components. *Note: Drawings are not to scale.*

## 4. Operation

### RondelO IV-set adapter – priming

1. Visually check packaging integrity.
2. Using aseptic technique, open the sterile packaging to access the RondelO IV-set adapter.
3. When applicable close the vented cap on the spike (1).
4. Verify the RondelO handle (7) is in a Flushing position.
5. Close the Pinch clamp (2) of the RondelO IV-set adapter.
6. Remove the protective covering from spike tip.

 The exposed spike tip should not come into contact with anything, otherwise it can become contaminated.

7. Spike the Infusion bag/container.
8. Twist off the spike port cap (8).
9. Of the full length IV-set to be connected, close the roller clamp and the vented cap on the spike if present, and remove the protective covering of the spike.
10. Spike the full length IV-set into the spike port (8) of the RondelO IV-set adapter.
11. Open the Pinch clamp (2) on the RondelO IV-set adapter.
12. Hold the Drip chamber tilted upside down and position it at the same level as the RondelO, open the Roller clamp of the IV-set and fill the Drip chamber as in (Fig 1.).
13. Continue to prime the rest of the IV-set.

 If all the air is not removed from the tubing, air may be forced into the patient's vascular system and cause an air embolism

14. Close the Roller clamp.
15. Connect the Primary IV-set to the patient's IV-access port.

 Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product

 All Luer lock connections should be finger tightened. Over tightening Luer lock connections can cause cracks and leaks to occur that could result in an air embolism and or exposure to biohazards.

### RondelO IV-set adapter – connect a secondary IV-set to the Drug position to be used of the RondelO

1. Verify the clamp on the secondary IV-set is closed before connecting.
2. Verify the RondelO handle (7) is not aligned with the Drug position.
3. If applicable, remove the protective cap (5) from the Luer lock connector at the Drug position.
4. Remove the cap from the Secondary IV-set and connect the Luer lock connector at the Drug position.
5. Repeat the procedure for each successive secondary IV-set.

### RondelO IV-set adapter – administer secondary infusions

1. Close the Pinch clamp (2) on the RondelO IV-set adapter above the RondelO.

2. Turn the RondelIO handle (7) to the Drug position and open the clamp on the secondary IV-set.
3. Use infusion pump to regulate flow rate to patient, alternatively by adjusting the roller clamp on the Primary IV-set.
4. Close the clamp on the secondary IV-set after the drug administration is completed.

### Flushing between and after drug administrations

1. Turn the RondelIO handle (7) to an adjacent Flushing position.
2. Open the Pinch clamp (2) above the RondelIO and flush.

### Disconnect from the RondelIO


1. Verify that the RondelIO handle (7) is not aligned with the Drug position to be disconnected.
2. Verify the clamp on the secondary IV-set is closed before disconnection.
3. Disconnect the secondary IV-set from the RondelIO.
4. If applicable, cap the Luer lock connector at the drug position of RondelIO.

### Disconnect the RondelIO IV-set adapter and the primary IV-set from the patient

1. Close the Roller clamp on the primary IV-set.
2. Disconnect it from the patient's IV access port.

## 5. RondelIO IV-set adapter, DC on top

**DEVICE DESCRIPTION:** The RondelIO IV-set adapter (DC on top) may consist of various components: (1) Vented spike with a protective covering integrated with Drip chamber, (2) Pinch clamp, (3) One-way valve, (4) RondelIO with 2, 4, or 6, Drug positions with either a (5) FLL cap, or (6) Needle free valve (NFV), (7) RondelIO Handle, (8) Spike port, Tubing, Transparent or UV-protected components. *Note: Drawings are not to scale.*

### RondelIO IV-set adapter (DC on top) – priming

1. Open the sterile packaging to access the RondelIO IV-set adapter.
2.  Product is sterile, do not use if package damaged or open.
3. When applicable, close the vented cap on the spike (1).
4. Verify the RondelIO handle (7) is in a Flushing position.
5. Close the Pinch clamp (2) of the RondelIO IV-set adapter.
6. Remove the protective covering from the spike tip.
7.  The exposed spike tip should not come into contact with anything, otherwise it can become contaminated.
8. Spike the infusion bag/ container.
9. Twist off the spike port cap (8).
10. Of the full length IV-set to be connected, close the roller clamp and the vented cap on the spike if present, and remove the protective covering.
11. Spike the full length IV-set into the adapter spike port (8) of the RondelIO IV-set adapter.
12. Open the Pinch clamp (2) of the RondelIO IV-set adapter.
13. Fill the Drip chamber (1) of the RondelIO IV-set adapter with the squeeze and release technique.
14. Hold the Drip chamber of the full length IV-set tilted upside down and position it at the same level as the RondelIO, open the Roller clamp of the IV-set and fill the Drip chamber as in (Fig 1.).

13. Continue to prime the rest of the IV-set.

If all air is not removed from the tubing, air may be forced into the patient's vascular system and cause an air embolism.

14. Close the Roller clamp.

15. Connect the Primary IV-set to the patient's IV-access port.

Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

All Luer lock connections should be finger tightened. Over tightening Luer lock connections can cause cracks and leaks to occur that could result in an air embolism and or exposure to biohazards.

### **RondelO IV- set adapter (DC on top) – connect a secondary IV-set to the Drug position to be used with the RondelO**

1. Verify the clamp on the secondary IV-set is closed before connecting.
2. Verify the RondelO handle (7) is not aligned with the Drug position.
3. If applicable, remove the protective cap (5) from the Luer Lock connector at the Drug position.
4. Remove the cap from the secondary IV-set and connect it to the Luer Lock connector at the Drug position.
5. Repeat the procedure for each successive secondary IV-set.

### **RondelO IV- set adapter (DC on top) – administer secondary infusions**

1. Utilise the head height differential practice or use the Pinch clamp on the RondelO IV-set adapter to regulate flow from primary or secondary infusion.
2. Turn the RondelO handle (7) to the Drug position and open the clamp fully on the secondary IV-set.
3. Use the infusion pump or adjust the roller clamp on the primary IV-set to regulate the flow rate.  
*Flushing fluid from the Primary bag will automatically take over after the Secondary IV-infusion is completed when using the head height differential practice.*

### **Flushing between and after drug administrations**

1. Turn the RondelO handle (7) to an adjacent Flushing position.
2. Verify or Open the Pinch clamp (2) above the RondelO and flush.

### **Disconnect from the RondelO**

1. Verify that the RondelO handle (7) is not aligned with the Drug position to be disconnected.
2. Verify the clamp on the secondary IV-set is closed before disconnection.
3. Disconnect the secondary IV-set from the RondelO.
4. If applicable, cap the Luer lock connector at the drug position of the RondelO.

### **Disconnect the RondelO IV-set adapter and the primary IV-set from the patient**

1. Close the Roller clamp on the primary IV-set.
2. Disconnect it from the patient's IV access port.

## **7. Storage and Disposal**

Keep dry

Handling mark on package – Fragile

Use by date

Keep away from sunlight

After having contact with the blood, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

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